

ZIANA INDUSTRIES

Diagnostic Instruments

CATALOGUE



Ziana Industries offers a wide range of laryngoscopes & Diagnostic Instruments. We have been in business since 1998 and providing our customers with best quality products and services. We welcome any inquiries/quotes that you may wish to raise please do not hesitate to contact us



ZIANA
INDUSTRIES

A B O U T U S

ZIANA SURGICAL is known for its quality for Surgical instruments and Diagnostic Instruments and is an operated business that has been dedicated to providing healthcare professionals with the products they need to perform their art since a long. Our high quality, precision crafted instruments have been an essential product line to our customers. We carry an exclusive range of Surgical Instruments and healthcare equipment's.

We are committed to providing a product line that guarantees quality and extraordinary value to the end user as we as a level of services, satisfaction, and cooperation that has become a hallmark of ZIANA INDUSTRIES.

We have gained a reputation for offering unbeatable prices on branded items as well as outstanding customers services on all levels. Every single instrument is manufactured with great care. We are working hard for maintain higher standard. Advanced techniques are used in manufacturing and finishing departments in order to made best products.

ZIANA Industries is offering you huge selection of our products and much more for a fast growing number of customers. We make it easy for you to browse our products, place order and benefit from our quality services.

S1

Section 01

1. Macintosh Blades – Standard
2. Miller Blades – Standard
3. Macintosh Blades – FIBEROPTIC
4. Macintosh Blades – FIBEROPTIC Fix Tube
5. Miller Blades - FIBEROPTIC
6. Miller Blades – FIBEROPTIC Fix Tube
7. MEGA BLADE - FIBEROPTIC

S2

Section 02

8. McCoy MEGA Blades – FIBEROPTIC
9. McCoy Macintosh Blades – Standards
10. McCoy Macintosh Blades – FIBEROPTIC
11. Miller Optima CLX Blades – FIBEROPTIC
12. Other Blades
16. Disposable Blades – Standard
17. Disposable Blades - FIBEROPTIC

S3

Section 03

18. Standard and F.O. Handles
20. Standard and F.O. Handles – New Pattern
22. Standard and F.O Handles - Rechargeable

S4

Section 04

24. Fiberoptic Laryngoscope Sets
25. Standard Laryngoscope

S5

Section 05

26. Otoscope & Ophthalmoscope Reusable
27. Otoscope & Ophthalmoscope Sets
30. Dermatoscope & ENT Set
31. Diagnostic Sets
32. Mini Otoscope
33. Mini Ophthalmoscope
34. Pocket Style Ophthalmoscope
35. Pocket Style Dermatoscope

S6

Section 06

36. Spare Accessories Diagnostic Set
38. Bulbs & Spatulas

S7

Section 07

39. Diagnostic Hammers

S8

Section 08

40. Tuning Forks

S9

Section 09

41. Rigid Laryngoscopes & Microsurgery

S10

Section 10

42. Magill Forceps



S11

Section 11

43. Case & Packing

CONTENT PAGE

ZIANA INDUSTRIES



Laryngoscopes

Conventional Blades
Fiberoptic Blades
American Profile Blades
Disposable Blades
Laryngoscope Handles
Laryngoscope Sets

ABOUT LARYNGOSCOPES

What is it?

A Laryngoscope is an instrument used to examine the larynx. The classic laryngoscope is available with a variety of blade attachments and with/without a Fiber Optic light source. The trend regarding blades has/is moving toward the disposable version.

Function

Laryngoscope blades are used as the primary tools for examination of the interior of the larynx and for placement of an endotracheal tube.

Main Feature

The Miller Laryngoscope Blade the most popular of the straight types. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-intubate patients.

The Macintosh Laryngoscope Blade the Macintosh is the predominate model of the curved blade types. During laryngoscopy, the tip or beak of the blade is compressed into the angle formed by the base of the tongue and the epiglottis, indirectly raising the epiglottis. Some sizes are available flangeless. Also available is a die-stamped or one-piece construction blade with a slightly expended flange.



FIBER OPTIC

Reusable Blade



- All blades in fiberoptic system crafted from finest non-magnet stainless steel minimize the risk of contamination.
- Cool and brighter light with no chance of tissue burning during long intubation procedures.
- Superior quality of glass fiber optic gives maximum bright illumination to view the larynx of human body.
- Removeable light guide can be autoclave separately, which makes long life of blade.
- Autoclave able at 134°C 5 minute, for 4000 cycles at maximum light output.
- All blades and handles are made to comply with ISO 7376, allow compatibility with existing blades and handles.
- Complete range of handles also equipped with LED lamps.
- LED technology has longer additional life time with 50,000hrs and continue 2hrs operating time.
- Brighter than Xenon lamp with 12,000 – 50,000 LUX transmission light at the tip of the blade.

STANDARD CONVENTIONAL Reusable Blade



We have accommodated the whole range of conventional lamp blades into LED choice also, which protects the patient from accidental burns during lengthy procedures.

- Standard Vacuum Lamp:
- Reliable vacuum standard lamp of 2.5V for better illumination.
- Conventional color temperature of 3000K @ 2000 LUX.
- 2hrs operating time with 20hrs life of lamp.

LED Lamp:

- LED lamp gives continuous cool and bright illumination.
- 2hrs operating time.
- Brighter than standard lamp with 12,000 LUX and also 85,000 LUX light beam.
- Available with 6500k color temperature.
- Longer additional life time with 50,000hrs.
- Whole range of blades complies with ISO 7376 standard.
- Easily replace lamp at the front of the blade.
- Autoclaved @ 134°C for 5 minutes.
- Crafted from finest stainless steel AISI 303 & 304.
- Round and smooth edges reduce the risk of injuries.

Macintosh Blades

Conventional Illumination

Specifications:


- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Available with LED and Standard Bulb Illumination.
- Innovation design for simple incubation and good visibility.
- Autoclave Able.
- Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux.
- These Blades Need No Maintenance –Ever But If Problem With Bulb Then It Can Be Replace.
- Quick-to-remove Lamp At The Front Of the Blades.
- The Rounded Edges of the Blades Reduce The Risk of Injury.
- **Require Conventional Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-C-0786	Premature	00	68 MM
ZI-C-0787	Neonate	0	80 MM
ZI-C-0788	Infant	1	95 MM
ZI-C-0789	Child	2	105 MM
ZI-C-0790	Medium	3	135 MM
ZI-C-0791	Large Adult	4	155 MM
ZI-C-0792	X-Large	5	175 MM



Section 01



Miller Blades


Conventional Illumination

Specifications:

- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Available with LED and Standard Bulb Illumination.
- Innovation design of the miller blade provides good visibility of the epiglottis and vocal cord.
- Autoclave Able.
- Bulbs 2.5v Halogen Also 3.5v Led With 20000 Lux.
- These Blades Need No Maintenance –Ever But If Problem With Bulb Then It Can Be Replace.
- Quick-to-remove Lamp At The Front Of the Blades.
- The Rounded Edges of the Blades Reduce The Risk of Injury.
- **Require Conventional Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-C-0793	Premature	00	65 MM
ZI-C-0794	Neonate	0	78 MM
ZI-C-0795	Infant	1	102 MM
ZI-C-0796	Child	2	154 MM
ZI-C-0797	Medium	3	195 MM
ZI-C-0798	Large Adult	4	205 MM



Section 01

Macintosh Blades



Fiber Optic Illumination

Specifications:

- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Autoclave Able.
- **Available In 3mm, 4mm, 5mm Fiber Optic Illumination.**
- Inserted Fiber Optic Technology For Advanced Performed
- Individual Micro-fibers For Improved Light Transmission and Longer Life.
- Smooth Shape Avoids Sharp Edges and Cavities For Easy Cleaning and Sterilizing.
- Extra-Strong **Removeable Fiber Bundles** To Ensure The Best Possible Light Intensity and Durability.
- The Clear Design to the Macintosh F.O. Blades simplifies Cleaning and Autoclaving.
- Fiber Optic Pipe is removable it can be remove or replace.
- **Require Fiber Optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0799	Premature	00	68 MM
ZI-F-0800	Neonate	0	80 MM
ZI-F-0801	Infant	1	95 MM
ZI-F-0802	Child	2	105 MM
ZI-F-0803	Medium	3	135 MM
ZI-F-0804	Large Adult	4	155 MM
ZI-F-0805	X-Large	5	175 MM



Section 01

Macintosh Blades


Fiber Optic Illumination

Specifications:

- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Autoclave Able.
- **Available In 3mm, 4mm, 5mm Fiber Optic Illumination.**
- Inserted Fiber Optic Technology For Advanced Performed
- Individual Micro-fibers For Improved Light Transmission and Longer Life.
- Smooth Shape Avoids Sharp Edges and Cavities For Easy Cleaning and Sterilizing.
- Extra-Strong **Fixed Fiber Bundles** To Ensure The Best Possible Light Intensity and Durability.
- The Clear Design to the Macintosh F.O. Blades simplifies Cleaning and Autoclaving.
- **Fix Tube** – Fiber Optic Pipe is Non-removeable it cannot be Remove.
- **Require Fiber Optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0806	Premature	00	68 MM
ZI-F-0807	Neonate	0	80 MM
ZI-F-0808	Infant	1	95 MM
ZI-F-0809	Child	2	105 MM
ZI-F-0810	Medium	3	135 MM
ZI-F-0811	Large Adult	4	155 MM
ZI-F-0812	X-Large	5	175 MM



Section 01

Miller Blade

Fiber Optic Illumination

Specifications:

- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Autoclave Able.
- **Available In 3mm, 4mm, 5mm Fiber Optic Illumination.**
- Fiber Optic Pipe is Removeable It can be Remove or Replace.
- Inserted Fiber Optic Technology Light Transmission and Longer life.
- Smooth Shape Avoids Sharp Edges and Cavities for Easy Cleaning and Sterilizing.
- The Flat Design of the F.O. Miller Blade Improves Visibility Considerably and Simplifies Incubation.
- This Blades Reduces The Risk of Trauma to the Patient and Provides Excellent Visualization.
- **Require Fiber Optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0813	Premature	00	65 MM
ZI-F-0814	Neonate	0	78 MM
ZI-F-0815	Infant	1	102 MM
ZI-F-0816	Child	2	154 MM
ZI-F-0817	Medium	3	195 MM
ZI-F-0818	Large Adult	4	205 MM



Section 01



Miller Blade

Fiber Optic Illumination

Specifications:

- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Autoclave Able.
- Available In 3mm, 4mm, 5mm Fiber Optic Illumination.**
- Fix Tube** - Fiber Optic Pipe is Non-removeable It cannot be Remove or Replace.
- Inserted Fiber Optic Technology for Advanced Performance.
- Individual Micro-fibers For Improved Light Transmission and Longer Life.
- Smooth Shape Avoids Sharp Edges and Cavities For Easy Cleaning and Sterilizing.
- The Flat Design of the F.O. Miller Blades Improves Visibility Considerably and simplifies incubation.
- These Blades Reduces the Risk of Trauma to the Patient and Provides Excellent Visualization.
- These Blades Need No Maintenance – Ever.
- There are no Connecting Screws to Get Lost, No External Light Pathways to get wrong
- Require Fiber Optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0819	Premature	00	65 MM
ZI-F-0820	Neonate	0	78 MM
ZI-F-0821	Infant	1	102 MM
ZI-F-0822	Child	2	154 MM
ZI-F-0823	Medium	3	195 MM
ZI-F-0824	Large Adult	4	205 MM



Section 01



Mega Blade

Fiber Optic Illumination

Specifications:

- **Premium Quality Blades.**
- Manufactured From AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard
- Autoclave Able
- Inserted Fiber Optic Technology for Advanced Performance.
- Up To 6500 Individual Micro-fibers For Improved Light Transmission and Longer Life.
- These Blades Need no Maintenance – Ever.
- 5 Years Warranty.
- Built With An integral, extra-large Fiber Optic Bundle Providing Coller, White Light.
- Smooth Shape Avoids sharp Edges and Cavities for Easy Cleaning and Sterilizing.
- There Are no connecting Screws to get lost, No External Light Pathways to get wrong.
- Both ends of the Blades are protected for life.
- **Require Fiber Optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0825	Neonate	0	80 MM
ZI-F-0826	Infant	1	92 MM
ZI-F-0827	Child	2	115 MM
ZI-F-0828	Medium	3	135 MM
ZI-F-0829	Large Adult	4	155 MM
ZI-F-0830	X-Large	5	175 MM



McCOY - Mega Blade

Fiber Optic Illumination

Specifications:

- The McCoy is a Macintosh Laryngoscope Blade With A Tip Which is adjustable through 70°. A lever Controls the Tip Angled During incubation to lift the Epiglottis, Giving a Clear View of the vocal Cords. Contact with the Upper Row of the Teeth and Unnecessary Pressure on the Tongue and soft Tissue are avoided.
- Manufactured From AISI 304 stainless steel.
- All Blades Compatible with ISO 7376 standard.
- Autoclave Able.
- Inserted Fiber Option Technology for Advanced Performance.
- Up To 65000 individual Micro-fibers for improved Light Transmission and Longer Life.
- These Blades Need no Maintenance – Ever.
- Built with An integral, Extra-Large Fiber Optic Bundle Providing Cooler, White Light.
- Smooth Shape Avoids sharp edges and cavities for easy cleaning and sterilizing.
- Integral Design Enables Superior Cleaning by Elimination of Cracks and Cavities.
- Both Ends of the Blades are protected for life.
- **Require Fiber Optic Handle.**



Section 02

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0831	Child	2	102 MM
ZI-F-0832	Medium	3	130 MM
ZI-F-0833	Large Adult	4	155 MM



McCOY - Blade

Standard Illumination

Specifications:


- The Mccoy Is a Macintosh Laryngoscope Blade With A tip which is adjustable Through 70°. A lever controls the tip angle during incubation of the lift the Epiglottis, Giving a Clear view of the vocal cords. Contact with the upper row of the teeth and unnecessary Pressure on the Tongue and soft Tissue are avoided.
- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Available with Led and standard Bulb Illumination.
- Innovation Design for simple incubation and good visibility.
- Autoclave Able.
- Bulbs 2.5v Halogen also 3.5v Led with 2000 Lux.
- These Blades need no maintenance – ever but if problem with bulb then it can be replace.
- Quick-to-remove Lamp at the front of the blade.
- The Rounded edges of the Blades reduce the risk of injury.
- Require Conventional Handle.



Section 02

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-C-0834	Child	2	102 MM
ZI-C-0835	Medium	3	130 MM
ZI-C-0836	Large Adult	4	155 MM



McCOY - Blade

Fiber Optic Illumination

Specifications:

- The McCoy Is a Macintosh Laryngoscope Blade With A tip which is adjustable Through 70°. A lever controls the tip angle during incubation of the lift the Epiglottis, Giving a Clear view of the vocal cords. Contact with the upper row of the teeth and unnecessary Pressure on the Tongue and soft Tissue are avoided.
- Manufactured from AISI 304 Stainless Steel.
- All Blades Compatible with ISO 7376 Standard.
- Autoclave Able
- Available in 3mm, 4mm, 5mm Fiber Optic Illumination
- Inserted Fiber Optic Technology for advanced performance.
- Individual Micro-fibers for Improved Light Transmission and Longer Life.
- Smooth Shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- Extra-strong Fiber Bundles to Ensure the Best Possible Light Intensity and Durability.
- The clear Design of the Macintosh F.O. Blades Simplifies cleaning and autoclaving.
- **Fix-Tube** – Fiber Optic pipe is non-removeable it can be remove.
- Also available in standard fiber optic Illumination (removeable pipe)
- **Require Fiber Optic Handle.**



Section 02

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0837	Child	2	102 MM
ZI-F-0838	Medium	3	130 MM
ZI-F-0839	Large Adult	4	155 MM



Miller Optima


Fiber Optic Illumination

Specifications:

- **Highly polished 4mm replaceable, glass fiber core for focused illumination. It has a one piece snap in/out fiber core for easy cleaning and replacement. No tools required.**
- Atraumatic distal end.
- Fully autoclave able metal blade with an anti-glare satin finish. Covered by a 5 years guarantee or 4,000 autoclave, whichever is sooner.
- A green, fiber optic ISO 7376 fitting metal block.
- 3 metal ball bearing for secure fitting.
- Clear product information etched on the blade.
- Available in 3mm, 4mm, 5mm Fiber Optic Illumination.
- Inserted Fiber Optic Technology for Advanced Performance.
- Individual Micro-fibers for improved Light Transmission and Longer Life.
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- Extra-strong Fiber Bundles to ensure the best possible light intensity and durability.
- Also available in standard illumination.
- **Require fiber optic Handle.**

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0840	Amature	000	60 MM
ZI-F-0841	Premature	00	66.5 MM
ZI-F-0842	Neonate	0	77 MM
ZI-F-0843	Infant	1	101 MM
ZI-F-0844	Child	2	153 MM
ZI-F-0845	Medium	3	197 MM
ZI-F-0846	Large Adult	4	205 MM



Section 02

American Profile - Blade

Macintosh Blades - Conventional Blades

Code	Blades	B/No	Overall Length
ZI-C-0847	Neonate	0	80 MM
ZI-C-0848	Infant	1	90 MM
ZI-C-0849	Child	2	108 MM
ZI-C-0850	Medium	3	135 MM
ZI-C-0851	Large Adult	4	160 MM



Section 02

Miller Blades – Conventional Blades

Code	Blades	B/No	Overall Length
ZI-C-0852	Premature	00	65 MM
ZI-C-0853	Neonate	0	75 MM
ZI-C-0854	Infant	1	102 MM
ZI-C-0855	Child	2	155 MM
ZI-C-0856	Medium	3	195 MM
ZI-C-0857	Large Adult	4	205 MM



American Profile - Blade

Philips Blades - Conventional Blades

Code	Blades	B/No	Overall Length
ZI-C-0858	Infant	1	113 MM
ZI-C-0859	Adult	2	166 MM



Wisconsin Blades – Conventional Blades

Code	Blades	B/No	Overall Length
ZI-C-0860	Child	2	132 MM
ZI-C-0861	Medium	3	160 MM
ZI-C-0862	Large Adult	4	200 MM



Polio Blades, English Profile – Conventional Blades

Code	Blades	B/No	Overall Length
ZI-C-0863	Child	2	105 MM
ZI-C-0864	Medium	3	130 MM
ZI-C-0865	Large Adult	4	155 MM



American Profile - Blade

Macintosh Blades – Fiber Optic Blades

Code	Blades	B/No	Overall Length
ZI-F-0866	Premature	00	68 MM
ZI-F-0867	Neonate	0	80 MM
ZI-F-0868	Infant	1	90 MM
ZI-F-0869	Child	2	108 MM
ZI-F-0870	Medium	3	130 MM
ZI-F-0871	Large Adult	4	160 MM



Section 02

Miller Blades – Fiber Optic Blades

Code	Blades	B/No	Overall Length
ZI-F-0872	Premature	00	65 MM
ZI-F-0873	Neonate	0	75 MM
ZI-F-0874	Infant	1	102 MM
ZI-F-0875	Child	2	155 MM
ZI-F-0876	Medium	3	195 MM
ZI-F-0877	Large Adult	4	205 MM



American Profile - Blade

Philips Blades – Fiber Optic Blades

Code	Blades	B/No	Overall Length
ZI-F-0878	Infant	1	113 MM
ZI-F-0879	Adult	2	166 MM



Wisconsin Blades – Fiber Optic Blades

Code	Blades	B/No	Overall Length
ZI-F-0880	Child	2	132 MM
ZI-F-0881	Medium	3	160 MM
ZI-F-0882	Large Adult	4	200 MM



Polio Blades, English Profile – Fiber Optic Blades

Code	Blades	B/No	Overall Length
ZI-F-0883	Child	2	105 MM
ZI-F-0884	Medium	3	130 MM
ZI-F-0885	Large Adult	4	155 MM



Disposable Blades

Standard Illumination

Specifications:

- Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- Packed individually in a clear pouch making it easy to identify blade size and type.
- Reduced risk of cross-contamination.
- Stainless steel body.
- Macintosh Sizes 0-4
- Millar sizes 00-4
- Interchangeable with all standards handle.



Size Measurement:

Code	Blades	B/No	Overall Length
ZI-C-0886	Premature	00	68 MM
ZI-C-0887	Neonate	0	80 MM
ZI-C-0888	Infant	1	95 MM
ZI-C-0889	Child	2	105 MM
ZI-C-0890	Medium	3	135 MM
ZI-C-0891	Large Adult	4	155 MM

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-C-0892	Neonate	0	80 MM
ZI-C-0893	Infant	1	95 MM
ZI-C-0894	Child	2	105 MM
ZI-C-0895	Medium	3	135 MM
ZI-C-0896	Large Adult	4	155 MM



Disposable Blades

Standard Illumination

Specifications:

- Single-use, sterile disposable blades help eliminate the time, reprocessing steps and costs related to disinfecting reusable blades.
- Packed individually in a clear pouch making it easy to identify blade size and type.
- Reduced risk of cross-contamination.
- Stainless steel body.
- Macintosh Sizes 0-4
- Millar sizes 00-4
- Interchangeable with all fiber optic handles.



Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0897	Premature	00	68 MM
ZI-F-0898	Neonate	0	80 MM
ZI-F-0899	Infant	1	95 MM
ZI-F-0900	Child	2	105 MM
ZI-F-0901	Medium	3	135 MM
ZI-F-0902	Large Adult	4	155 MM

Size Measurement:

Code	Blades	B/No	Overall Length
ZI-F-0903	Neonate	0	80 MM
ZI-F-0904	Infant	1	95 MM
ZI-F-0905	Child	2	105 MM
ZI-F-0906	Medium	3	135 MM
ZI-F-0907	Large Adult	4	155 MM



Handles

Fiber Optic Illumination

Specifications:

- Chrome Plated Brass Fiber Optic Handle compatible with all "Green System" fiber optic blades.
- Require 2 AA size Batteries or 3.5 V Rechargeable Battery.
- Crafted from high quality Chrome Plated Brass.
- Autoclavable.

Replacement Bulbs

Bulb	Battery
Xenon	2.5v
LED	3.0v
Halogan	2.5v
Halogan	3.5v

Battery Handle Fiberoptic Xenon, Brass Chrome Plated

Code	Battery	Handle	Overall Length
ZI-F-0908	1 x AA Size	Infant	122 MM
ZI-F-0909	1 x C Size	Stubby	125 MM
ZI-F-0910	2 x AA Size	Pediatric	160 MM
ZI-F-0911	2 x C Size	Medium	160 MM
ZI-F-0912	2 x D Size	Large	175 MM

Battery Handle Fiberoptic LED, Brass Chrome Plated

Code	Battery	Handle	Overall Length
ZI-F-0913	1 x AA Size	Infant	122 MM
ZI-F-0914	1 x C Size	Stubby	125 MM
ZI-F-0915	2 x AA Size	Pediatric	160 MM
ZI-F-0916	2 x C Size	Medium	160 MM
ZI-F-0917	2 x D Size	Large	175 MM



Handles

Conventional Illumination

Specifications:

- Chrome Plated Brass Conventional Handle compatible with all Conventional blades.
- Require 2 AA size Batteries or 3.5 V Rechargeable Battery.
- Crafted from high quality Chrome Plated Brass.
- Autoclavable.

Replacement Bulbs

Bulb	Battery
Xenon	2.5v
LED	3.0v
Halogan	2.5v
Halogan	3.5v

Battery Handle Conventional, Brass Chrome Plated

Code	Battery	Handle	Overall Length
ZI-C-0918	1 x AA Size	Infant	110 MM
ZI-C-0919	1 x C Size	Stubby	115 MM
ZI-C-0920	2 x AA Size	Pediatric	145 MM
ZI-C-0921	2 x C Size	Medium	150 MM
ZI-C-0922	2 x D Size	Large	165 MM





Specifications:

- Chrome Plated Brass Half Plain and Double Round Cutting in Middle.
- Available standard and fiberoptic.
- Require 2 AA Size Batteries or 3.5 V Rechargeable battery.
- Crafted from High Quality Chrome plated brass.
- Autoclavable

Fiberoptic & Standard Battery Handle

Brass Chrome Plated

Code	Handle	Overall Length
ZI-F-0923	Infant	122 MM
ZI-F-0924	Stubby	125 MM
ZI-F-0925	Pediatric	160 MM
ZI-F-0926	Medium	160 MM
ZI-F-0927	Large	175 MM

Specifications:

- Round Handle with Circular Pattern From Top to Bottom.
- Available Standard and Fiberoptic.
- Require 2 AA Size Batteries or 3.5 V Rechargeable Battery.
- Available in Brass and Steel.
- Autoclave Able

Fiberoptic & Standard Battery Handle

Brass Chrome Plated

Code	Handle	Overall Length
ZI-F-0928	Infant	122 MM
ZI-F-0929	Stubby	125 MM
ZI-F-0930	Pediatric	160 MM
ZI-F-0931	Medium	160 MM
ZI-F-0932	Large	175 MM





Specifications:

- High Quality Battery Handle Round Circle on Bottom Top is Quality Symbol.
- Available Standard and Fiberoptic.
- Available in Brass and Steel.
- Autoclave able.

**Fiberoptic & Standard Battery Handle
Brass Chrome Plated**

Code	Handle	Overall Length
ZI-F-0933	Infant	122 MM
ZI-F-0934	Stubby	125 MM
ZI-F-0935	Pediatric	160 MM
ZI-F-0936	Medium	160 MM
ZI-F-0937	Large	175 MM

Specifications:

- High Quality Battery Handle with Half Learning Style.
- Available Standard and Fiberoptic.
- Available in Brass and Steel.
- Autoclave Able

**Fiberoptic & Standard Battery Handle
Brass Chrome Plated**

Code	Handle	Overall Length
ZI-F-0938	Infant	122 MM
ZI-F-0939	Stubby	125 MM
ZI-F-0940	Pediatric	160 MM
ZI-F-0941	Medium	160 MM
ZI-F-0942	Large	175 MM





Section 03

Handle

Fiber Optic Illumination

Code	Handle	Overall Length
ZI-F-0943	Pediatric	160 MM
ZI-F-0944	Medium	150 MM
ZI-F-0945	Medium with DC	150 MM

Handle

Conventional Illumination

Code	Handle	Overall Length
ZI-C-0946	Pediatric	145 MM
ZI-C-0948	Medium	150 MM
ZI-C-0949	Medium with DC	150 MM



Handle
Otoscope & Ophthalmoscope

Code	Handle	Overall Length
ZI-F-0950	Aluminum Handle	145 MM
ZI-F-0951	Medium	150 MM

Specifications:

- Chrome Plated Brass Fiber Optic Handle Compatible with all Fiber Optic Blades.
- 3.5 V Rechargeable Battery.
- Crafted from High Quality Chrome Plated Brass.
- Autoclave Able

Laryngoscope Sets

Fiber Optic Illumination



FO Macintosh BLADE SETS

Standard Sets

Code	Contents
ZI-F-0952	Set of 3 Blades Size 0, 1, 2 and Adult Handle
ZI-F-0953	Set of 3 Blades Size 1, 2, 3 and Adult Handle
ZI-F-0954	Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle
ZI-F-0955	Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle



FO MILLER BLADE SETS

Standard Sets

Code	Contents
ZI-F-0956	Set of 3 Blades Size 00, 0, 1 and Child Handle
ZI-F-0957	Set of 3 Blades Size 0, 1, 2 and Child Handle
ZI-F-0958	Set of 4 Blades Size 00, 0, 1, 2 and Child Handle
ZI-F-0959	Set of 4 Blades Size 0, 1, 2, 3 and Child Handle



FO MCCOY BLADE SETS

Standard Sets

Code	Contents
ZI-F-0960	Set of 2 Blades Size 2 and Adult Handle
ZI-F-0961	Set of 3 Blades Size 3 and Adult Handle
ZI-F-0962	Set of 4 Blades Size 2, 3 and Adult Handle
ZI-F-0963	Set of 4 Blades Size 2, 3, 4 and Child Handle

- Set Available in Hard Case or Soft Leather Pouch.
- Handle Option Available.
- Custom Set Available.
- For Sets in pouch **Add -P** after Code.

Laryngoscope Sets

Standard Illumination

MACINTOSH CONVENTIONAL SETS

Standard Sets

Code	Contents
ZI-C-0964	Set of 3 Blades Size 0, 1, 2 and Adult Handle
ZI-C-0965	Set of 3 Blades Size 1, 2, 3 and Adult Handle
ZI-C-0966	Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle
ZI-C-0967	Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle



MILLER CONVENTIONAL SETS

Standard Sets

Code	Contents
ZI-C-0968	Set of 3 Blades Size 00, 0, 1 and Child Handle
ZI-C-0969	Set of 3 Blades Size 0, 1, 2 and Child Handle
ZI-C-0970	Set of 4 Blades Size 00, 0, 1, 2 and Child Handle
ZI-C-0971	Set of 4 Blades Size 0, 1, 2, 3 and Child Handle



MCCOY CONVENTIONAL SETS


Standard Sets

Code	Contents
ZI-C-0972	Set of 2 Blades Size 2 and Adult Handle
ZI-C-0973	Set of 3 Blades Size 3 and Adult Handle
ZI-C-0974	Set of 4 Blades Size 2, 3 and Adult Handle
ZI-C-0975	Set of 4 Blades Size 2, 3, 4 and Child Handle



- Set Available in Hard Case or Soft Leather Pouch.
- Handle Option Available.
- Custom Set Available.
- For Sets in pouch **Add -P** after Code.

ZIANA INDUSTRIES



ENT DEVICES

Diagnostic Sets

Ophthalmoscope

Otoscope

Dermatoscope

Mini Otoscope & Ophthalmoscope

ABOUT OTOSCOPE

What is it?

An otoscope is a tool which shines a beam of light to help visualize and examine the condition of the ear canal and eardrum. Examining the ear can reveal the cause of symptoms such as an earache, the ear feeling full, or hearing loss.

Function

An otoscope or auriscope is a medical device which is used to look into the ears. Health care providers use otoscopes to screen for illness during regular check-ups and also to investigate ear symptoms. An otoscope potentially gives a view of the ear canal and tympanic membrane or eardrum.

Otoscope Features



Advanced LED
Illumination



Optical Quality 3X
Magnification Lens



Spare Specula



360 Rotating Lens

ABOUT OPHTHALMOSCOPE


What is it?

The purpose of ophthalmoscopy is to examine the back of the eye and make sure that it is healthy. Ophthalmoscopy is a routine part of a comprehensive eye exam. Healthcare providers also may perform ophthalmoscopy during a routine physical.

Function

It is used to detect and evaluate symptoms of retinal detachment or eye diseases such as glaucoma. Ophthalmoscopy may also be done if you have signs or symptoms of high blood pressure, diabetes, or other diseases that affect the blood vessels.

Main Feature



ABOUT OPHTHALMOSCOPE

What is it?

A dermatoscope can assess structures to the depth of reticular dermis, and record images for future comparison. The basic principle of dermoscopy is transillumination of a lesion in order to study it with high magnification to visualize subtle features.

Function

Dermoscopy is performed with a handheld instrument called a dermatoscope. The procedure allows for the visualization of subsurface skin structures in the epidermis, at the dermoepidermal junction, and in the papillary dermis; these structures are usually not visible to the naked eye [2-4].



Otoscopes & Ophthalmoscopes

Reusable



ZI-O-0976
Brass Oscope Fiber Optic
Brass Oscope Conventional



ZI-O-0977
ENT Complete Set
Pouch Packing



ZI-O-0978
Brass
Ophthalmoscope

Specifications:

- Manufactured from Brass Chrome Plated.
- 3x Magnifying glass.
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles require 2 x C size Batteries.
- Hard Case and Soft Pouch Available.
- Fiber Optic Light Source.

Fitting



Pin Fitting



Bayonet Fitting

Otoscopes & Ophthalmoscopes

Standard Sets

BASIC OTOSCOPE SET (Bayonet Fitting)

Complete with Otoscope Battery Handle, 1 set of 5 Reusable specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm in a Pouch

Code	Contents
ZI-O-0979	Complete Set with Spare Lamp
ZI-O-0980	Complete set without Spare Lamp
ZI-O-0981	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-0982	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-0983	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



Section 05

BASIC OPHTHALMOSCOPE SET

Complete with Ophthalmoscope Battery Handle in a Pouch

Code	Contents
ZI-O-0984	Complete Set with Spare Lamp
ZI-O-0985	Complete set without Spare Lamp
ZI-O-0986	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-0987	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-0988	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



Otoscopes & Ophthalmoscopes

Standard Sets

STANDARD OTOSCOPE SET (Pin Fitting)

Complete with Otoscope Battery Handle, 1 set of 5 Reusable specula Dia 2.5, 3.5, 4.5, 5.5, 9.0mm in a Pouch

Code	Contents
ZI-O-0989	Complete Set with Spare Lamp
ZI-O-0990	Complete set without Spare Lamp
ZI-O-0991	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-0992	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-0993	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



Section 05

BASIC NASAL OTOSCOPE SET (Pin Fitting)

Complete with Otoscope Battery Handle, 1 Nasal Speculum, In a Pouch also Available in Plastic Box

Code	Contents
ZI-O-0994	Complete Set with Spare Lamp
ZI-O-0995	Complete set without Spare Lamp
ZI-O-0996	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-0997	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-0998	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



Otosopes & Ophthalmoscopes

Standard Sets

BASIC OPHTHALMOSCOPE SET

Complete with Ophthalmoscope Battery Handle, In A Plastic Box

Code	Contents
ZI-O-0999	Complete Set with Spare Lamp
ZI-O-1000	Complete set without Spare Lamp
ZI-O-1001	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-1002	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-1003	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



STANDARD OTOSCOPE & OPHTHALMOSCOPE

Complete with Otoscope & Ophthalmoscope Battery Handle, In A Plastic Box

Code	Contents
ZI-O-1004	Complete Set with Spare Lamp
ZI-O-1005	Complete set without Spare Lamp
ZI-O-1006	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-1007	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-1008	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



STANDARD OTOSCOPE SET

Complete with Otoscope, Battery Handle, 1 Set of 3 Reusable Specula Dia 4.5, 5.5, 9.0mm In A Plastic Box

Code	Contents
ZI-O-1009	Complete Set with Spare Lamp
ZI-O-1010	Complete set without Spare Lamp
ZI-O-1011	Complete Set with Spare Lamp and Rechargeable Battery
ZI-O-1012	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-O-1013	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



BASIC DERMATOSCOPES SET

Complete with Dermatoscopes, Spare Mirror In A Plastic Box


Code	Contents
ZI-D-1014	Complete Set with Spare Lamp
ZI-D-1015	Complete set without Spare Lamp
ZI-D-1016	Complete Set with Spare Lamp and Rechargeable Battery
ZI-D-1017	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-D-1018	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



BASIC DERMATOSCOPES SET

Complete With Dermatoscopes, Spare Mirror, In A Pouch

Code	Contents
ZI-D-1019	Complete Set with Spare Lamp
ZI-D-1020	Complete set without Spare Lamp
ZI-D-1021	Complete Set with Spare Lamp and Rechargeable Battery
ZI-D-1022	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-D-1023	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



Section 05

STANDARD DIAGNOSTIC SET

Complete with Standard Otoscope, Standard Ophthalmoscope, Handle, Nasal Speculums, Laryngeal Mirror, Reusable Specula, Hard Case

Code	Contents
ZI-D-1024	Complete Set with Spare Lamp
ZI-D-1025	Complete set without Spare Lamp
ZI-D-1026	Complete Set with Spare Lamp and Rechargeable Battery
ZI-D-1027	Complete Set with Spare Lamp, Rechargeable Battery & Desk Charger
ZI-D-1028	Complete Set with Spare Lamp, Rechargeable Battery & Wall Charger



DIAGNOSTIC SET

Full Specification



UNIVERSAL DIAGNOSTIC SET

ZI-D-1029

- Handle (dry Cell or Rechargeable) X 1
- Otoscope
- Ophthalmoscope
- Laryngeal Mirrors X2
- Reusable Specula X5
- Tongue Blade Holder
- Metal Tongue Blade Depressors
- Nasal Speculum
- Ear Tip Connector
- Ear Tips
- Replacement Bulb X 2
- Hard Case X 1
- Also Available in Leather Pouch Just add **-P** after code.



STANDARD DIAGNOSTIC SET

ZI-D-1030

- Standard Otoscope X 1
- Standard Ophthalmoscope X 1
- Handle (dry Cell Or Rechargeable) X 1
- Nasal Speculum
- Laryngeal Mirrors 2x
- Reusable Specula X 5
- Replacement Bulb X 2
- Hard Case X 1
- Also Available in Leather Pouch Just add **-P** after code.

About Sets

Diagnostic Sets are the Most Commonly Used Tools Sets During a clinical repertories to perform a Comprehensive Physical Test. Diagnostic Sets Facilitate Physicians to Perform Quick Exams & Accurately.

More Details

- Manufactured from High Quality Material.
- Available in Halogen & LED Illumination
- Bulbs 2.5V
- Handles Requires 2 x C size batteries
- Hard Case and Soft Pouch Available



ZI-O-1031

Mini Oscope Fiber Optic
Complete Set of Disposable Cannulas



ZI-O-1032

Mini Oscope Fiber Optic
Complete Set of Disposable Cannulas

Specification:

- Manufactured from Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA Size Batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source



Color Options





ZI-O-1033



ZI-O-1034

Mini Ophthalmoscope

Parker mini ophthalmoscope with 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses – 25 to + 40 for excellent resolution (Red for – Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial clear identification for adopter settings.

Power Supply: 2 x AA Batteries (not included in set).

Packing: Hard Plastic Case Packing.



ZI-O-1035



ZI-O-1036

Section 05



ZI-O-1037



Multi-Action Tongue Depressor

All-Purpose Pocket-Sized Illumination that is easy to use and maintain. Multiple Applications. Tongue-Depressor Holder can be used with disposable depressor (Included) or with Standard Wooden Depressors. Bright Halogen Bulb with More Intensity Compared with Conventional Bulbs. Refined Plastic Compact Body for Long Lasting Use. Attachment clip with on/off Button. Set Includes: (1) Handle, (1) Head with Integrated Tongue-Depressor Holder (3) Disposable Depressors LED Illumination Probe. Pouch. Two “AA” Batteries Required.

Batteries not included.



The Ophthalmoscope is a tool used in medical to examine the interior of the eye. Studying the inner parts of the eye, such as the optic disk, retinal blood vessels, retina, choroid and macula, can help diagnose diseases of the eye.

Pocket Ophthalmoscope

ZI-O-1038

The Pocket size ophthalmoscope is a lightweight, compact diagnostic Instrument that delivers high performance in a slightly smaller package. Quality optics in mini-format is deal for rounds and house calls. Attachment Clip with integrated on/off switch. High-quality plastic handle with shockproof, sturdy, non-slip features, Enhanced LED technology for perfect illumination and a brilliant image. 4 different apertures. Comprehensive lens range from 20 D to -20 D.

Power Supply: 2 X AA Batteries (not included in set).
Packing: Hard plastic case packing.

Pocket Ophthalmoscope

ZI-O-1039

2.5v LED lamp light for true tissue color and consistent long lasting illumination 16 colors coded lenses -25 to +25 for excellent resolution (Red for – Black for +) 4 different apertures (small, large, blue, half moon focus). Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification of adopter setting.

Power Supply: 2 X AA Batteries (not included in set).
Packing: Hard Plastic case packing.





Pocket Dermatoscope

ZI-D-1040

Modern pocket Dermatoscope in a compact design with high-quality optics. 10X magnification large, high-resolution image. Choice of two skin-friendly and sterilizable contact plates: with scale from 0 – 10mm for exact measurement of pigmented skin lesions, without scale. Simple lamp change on the underside of the instruments head. Enhanced LED Technology for more light. Attachment clip with integrated On/off switch. High-quality 2 color handle.

Power Supply: 2 X AA Batteries (Not included in set).

Packing: Hard Plastic case packing.



Dermatoscope

ZI-D-1041

Bayonet locking system for easy engagement and disengagement from the handle. 2.5v LED illumination of the examination field. High Resolution achromatic lenses system with 10x magnification. Adjustable focusing Eyepiece for a large correction range. Lens without scale. Soft eyepiece protection.



Available in Polish and Sand Finish.

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard Plastic case packing.



ZIANA INDUSTRIES



ACCESSORIES

Ear Specula
Testing Tube
Nasal
Speculum
Bulbs
Hammers
Tunning Forks
Forceps
Packing Case

TONGUE DEPRESSOR

Code Description

ZI-T-1042 Tongue Depressor Thermoplastic

METAL TONGUE DEPRESSOR

Code Description

ZI-T-1043 Tongue Depressor Metal

OTOSCOPE HEAD

Code Description

ZI-O-1044 Without Ear Spatula

ZI-O-1045 With Ear Spatula

ZI-O-1046 Bayonet Fitting

ZI-O-1047 Screw Fitting

ZI-O-1048 With LED

ZI-O-1049 With Bulb

ADAPTER FOR DISPOSABLE T.D

Code Description

ZI-A-1050 Adapter for Disposable T.D

METAL ADAPTER FOR SPECULA

Code Description

ZI-O-1051 Metal Adapter for BY Fitting

ZI-O-1052 Metal Adapter for PN Fitting

INSUFFLATOR BULB & TUBE

Code Description

ZI-B-1053 Insufflator Bulb & Tube



CURVED ILLUMINATOR



Code	Description
ZI-I-1054	Bayonet Fitting
ZI-I-1055	Pin Fitting
ZI-I-1056	Reusable

OPHTHALMOSCOPE HEAD



Code	Description
ZI-O-1057	Bayonet Fitting
ZI-O-1058	Pin Fitting
ZI-O-1059	Reusable

NASAL SPECULUMS



Code	Description
ZI-N-1060	Nasal Speculum



SPREADABLE NASAL SPECULUMS

Code	Description
ZI-N-1061	Spreadable Nasal Speculum



LARYNGEAL MIRRORS

Code	Description
ZI-L-1062	Laryngeal Mirrors



TONGUE DEPRESSOR

Code	Description
ZI-T-1063	Tongue Depressor with Bulb
ZI-T-1064	Tongue Depressor with LED

- Made from Thermosplastic Material.
- Smooth & Shiny Surface Ensure Easy Cleaning.
- Bayonet Lock for Otoscope.

DISPOABLE SPECULA

Code	Description
ZI-S-1065	Disposable Specula Dia 2.5mm, Box of 1000
ZI-S-1066	Disposable Specula Dia 4.0mm, Box of 1000
ZI-S-1067	Disposable Specula Dia 2.5mm, Box of 250
ZI-S-1068	Disposable Specula Dia 4.0mm, Box of 250
ZI-S-1069	Dia 2.5mm & 4.0mm in Tube
ZI-S-1070	Clear Tube for Specula



- Made from Thermosplastic Material.
- Smooth & Shiny Surface Ensure Easy Cleaning.
- Autoclave able at 130°C.
- Metallic Bayonet Lock for Otoscope.
- Can be disinfected.

REUSABLE SPECULA

Code	Description
ZI-S-1071	Reusable Specula No. 1 Dia 2.5mm
ZI-S-1072	Reusable Specula No. 2 Dia 3.5mm
ZI-S-1073	Reusable Specula No. 3 Dia 4.5mm
ZI-S-1074	Reusable Specula No. 4 Dia 5.5mm
ZI-S-1075	Reusable Specula No. 5 Dia 9.5mm



- Made from Thermosplastic Material.
- Smooth & Shiny Surface Ensure Easy Cleaning.
- Autoclave able at 130°C.
- Metallic Bayonet Lock for Otoscope.
- Can be disinfected.

REUSABLE VETERINARY SPECULA

Code	Description
ZI-S-1076	Reusable Veterinary Specula No. 1 Dia 4mm
ZI-S-1077	Reusable Veterinary Specula No. 2 Dia 5mm
ZI-S-1078	Reusable Veterinary Specula No. 3 Dia 7mm



Section 06




For Cheap conventional Laryngoscopes Blades
LUX : 16,000

Code	Item
ZI-B-1079	Bulb W/Lining
ZI-B-1080	Bulb W/O Lining



For conventional Laryngoscopes Blades
LUX : 20,000

Code	Item
ZI-B-1081	Bulb W/Lining
ZI-B-1082	Bulb W/O Lining



For Fiber Optic Laryngoscopes Blades
LUX: 20,000 to 40,000

Code	Item
ZI-L-1083	LED 20,000L
ZI-L-1084	LED 40,000L



Best Quality LED with Almost 50,000 LUX

Code	Item
ZI-L-1085	LED 35,000L
ZI-L-1086	LED 50,000L



LED Capsule Fix with Almost 50,000 to 100,000 LUX

Code	Item
ZI-L-1087	LED 50,000L
ZI-L-1088	LED 100,000L



For Mini & Metal Body Ophthalmoscope

Code	Item
ZI-B-1089	Bulb
ZI-L-1090	LED



For 4 Aperture Ophthalmoscope

Code	Item
ZI-B-1091	Bulb
ZI-L-1092	LED



Xenon

Code	Item
ZI-X-1093	Xenon



Halogen

Code	Item
ZI-H-1094	Halogen 2.5V
ZI-H-1095	Halogen 3.5V



Standard Capsule with Replaceable Bulb & LED

Code	Item
ZI-B-1096	With Bulb
ZI-L-1097	With LED

TAYLOR HAMMER

Code	Description
ZI-H-1098	Taylor Hammer S.S, 17 cm
ZI-H-1099	Taylor Hammer S.S, 20 cm

- The most popular Hammer
- Stainless Steel Handle
- Latex Free
- Available in Two Size



TROMNER HAMMER

Code	Description
ZI-H-1100	Tromner Hammer, 175Gm, S.S, 24cm

- S.S Handle with 2 Sided Synthetic Rubber Head
- Latex Free
- Available in 24 cm



BERLINER HAMMER

Code	Description
ZI-H-1101	Berliner Hammer S.S, 17 cm

- Stainless Steel Handle also with plastic grip.
- Latex Free
- Available Size 17 cm



DEJERINE HAMMER

Code	Description
ZI-H-1102	Dejerine Hammer, S.S, 20cm
ZI-H-1103	Dejerine Hammer with Needle, S.S, 21cm

- S.S handle with 2 sided synthetic Rubber Head
- Latex Free
- Handle Conceals Screw-in-bush
- Needle contained within the head



BUCK HAMMER

Code	Description
ZI-H-1104	Buck Hammer, S.S, 15cm
ZI-H-1105	Buck Hammer, S.S, 18cm

- S.S Handle with 2 sided synthetic Rubber Head
- Latex Free
- Handle Conceals Screw-in-bush
- Needle Contained within the head



BABINSK HAMMER

Code	Description
ZI-H-1106	Babinski Hammer with Needle, S.S, 20cm

- Octagonal Stainless Steel Handle with Concealed Needle
- Rubber Head Can only used in one Direction
- Latex Free



BABINSK HAMMER

Code	Description
ZI-H-1107	Babinski Hammer with Needle, S.S, 20cm

- Octagonal Stainless Steel Handle with Concealed Needle
- Rubber Head Can only used in one Direction
- Latex Free



Tuning Fork

Full Specification



HARTMAN



WELDON ALLOY



WELDON ALLOY
With Weight



KONIG



LUCE



RYDEL SEIFFERT

Code	Hz	Code	Hz	Code	Hz	Code	Hz	Code	Hz	Code	Hz
ZI-T-1108	128	ZI-T-1114	128	ZI-T-1120	128	ZI-T-1122	128	ZI-T-1123	256	ZI-T-1125	64
ZI-T-1109	256	ZI-T-1115	256	ZI-T-1121	256			ZI-T-1124	512		
ZI-T-1110	512	ZI-T-1116	512								
ZI-T-1111	1024	ZI-T-1117	1024								
ZI-T-1112	2048	ZI-T-1118	2048								
ZI-T-1113	4096	ZI-T-1119	4096								

Section 08



KONIG

Specifications:


- Manufactured from Aluminum Also avail in steel.
- Pouch Packing Available for single piece.
- Set available in hard case or Plastic Box or Wooden Box
- Custom Set Available**
- All tuning forks available in 6 different color.

ABOUT TUNING FORKS

- The Fundamental Mode of Vibrations is the Mode Most Commonly Associated with Tuning forks**
- The two tines of the Fork Alternately Move Toward and away from Each Other, Each Bending Like A Cantilever Beam, Fixed at the Stem and Free at the other end. This is a symmetric mode, Since the two Tines are Mirror Images of Each other



Code	Hz
ZI-T-1126	128



Jackson Laryngoscope

With handle and light carrier, stain finish inside

Size 1 (12.0cm)
ZI-L-1127

Size 2 (17.1cm)
ZI-L-1128

Size 3 (18.0cm)
ZI-L-1129



Jackson Laryngoscope

With handle and light carrier, stain finish inside

Size 1 (12.0cm)
ZI-L-1130

Size 2 (17.1cm)
ZI-L-1131

Size 3 (18.0cm)
ZI-L-1132



KLEINSASSER LARYNGOSCOPE

Anterior commissure, 22 cm especially for long neck

Size 1 (13.0cm)
ZI-L-1133

Size 2 (15.0cm)
ZI-L-1134


Size 3 (17.0cm)
ZI-L-1135

Size 3 (22.0cm)
ZI-L-1136


Rigid Laryngoscope




Esophgoscope with Fiber Optic Light Carrier
ZI-E-1137 25cm




Esophgoscope
ZI-E-1138 35cm



Bronchoscope
ZI-B-1139 45cm



Negus Aspirating Bronchoscope
ZI-B-1140 52cm



Adolescent Bronchoscope
ZI-B-1141 25cm
ZI-B-1142 30cm
ZI-B-1143 40cm

Magill Forceps

Full Specification



ZI-F-1144
17 cm




ZI-F-1145
20 cm



ZI-F-1146
25 cm

About Forceps:

- ABent Blunt Instruments used to facilitate Nasotracheal Intubation.
- Angled Forceps used to Guide a Tracheal Tube into the Larynx or a Nasogastric Tube into the Esophagus Under Direct Vision. It is also used to place Pharyngeal Packs and remove foreign bodies.



Packing



Made from Thermoplastic.
For All Kinds of Mac, Miller & Diagnostic Sets.
Available in 6 Colors.




Made from Thermoplastic.
For All Kinds of Mac, Miller & Diagnostic Sets.
Available in 6 Colors.



Specially For Mini Otoscope & Ophthalmoscope. Also available for all kinds of handle.
Available in All Colors

Box Packing



ZIANA INDUSTRIES

ZIANA INDUSTRIES

Address: Islam Nagar Road Pacca Garah,
Sialkot-51310, Pakistan.

Ph # +92 52 4261495 / +92 52 4265791

Fax # +92 52 4264970

Mobile/Whatsapp # +92 321 9408040

Email: sales@zianaindustries.com / info@zianaindustries.com /
expziana@gmail.com

URL: www.zianaindustries.com



Diagnostic Instruments